Serum dextrose agar (SDA)

Détails

For 1 liter of SDA:

33 g of Difco Tryptose blood agar base	or	10 g of Tryptose 3 g of Beef extract 5 g of NaCl 15 g of Agar

Complete to 950 ml with water.
Autoclave.

Add 50 ml of sterile horse serum.
Add 50 ml g of sterile glucose 20%.

Reference:
Corbel MJ, Thomas EL. 1980.
The Brucella-phages: their properties, characterization and applications. Booklet 2266, Ministry of
Agriculture, Fisheries and Food, Pinner, Middlesex, UK, 106 p

Last modifications: 2011-02-24