

## Media recipe

## **Albimi Brucella Broth**

## Détails

For 1 liter of Albimi:

20 g of peptone.

1 g of glucose.

0.1 g of sodium bisulphite.

2 g of yeast extract.

5 g of NaCl.

Complete to 1 liter with water.

For Albimi-agar:

Add 15 g of agar per liter.

For Albimi-top agar:

Add 7.5 g of agar per liter.

Autoclave.

Corbel MJ, Thomas EL. 1980. The Brucella-phages: their properties, characterization and applications. Booklet 2266, Ministry of Agriculture, Fisheries and Food, Pinner, Middlesex, UK, 106 p.

Last modifications: 2011-01-24