

Media recipe

AMS

Détails

For 1 liter of GM17:

8 g of NB (nutrient broth).

0.7 g of KCl.

1.5 g of $\text{CaCl}_2 \cdot 2\text{H}_2\text{O}$.

12 g of $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$.

26 g of NaCl

Complete to 1 liter with water.

For AMS-agar:

Add 15 g of agar per liter.

For AMS-top agar:

Add 7.5 g of agar per liter.

Autoclave.

Last modifications: 2011-01-16