

GAM

Détails

For 1 liter of GAM:

- 10 g of peptone.
- 3 g of soybean peptone.
- 10 g of proteose peptone.
- 13 g of digest serum powder.
- 5 g of yeast extract.
- 2,2 g of meat extract.
- 1,2 g of liver extract.
- 3 g of glucose.
- 2.5 g of dipotassium phosphate.
- 3 g of NaCl.
- 5 g of soluble starch.
- 0.3 g of L-cysteine-HCl.
- 0.3 g of sodium thioglycolate.
- Complete to 1 liter with water.

For GAM-agar:

Add 12g of agar per liter.

For GAM-top agar:

Add 7.5 g of agar per liter.

Autoclave.

Last modifications: 2011-01-24