

SWC

For 1 liter of SWC:

5 g of Peptone.

0.5 g of Yeast extract.

3 ml of Glycerol.

750 ml of Sea water.

Complete to 1 liter with water.

Autoclave.

Sea water

For 1 liter of Sea water:

23.376 g of NaCl (400 mM).

24.647 g of $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$ (100 mM).

1.49 g of KCl (20 mM).

2.94 g of $\text{CaCl}_2 \cdot 2\text{H}_2\text{O}$ (20 mM).

Complete to 1 liter with water.

Autoclave.

Last modifications: 2011-02-24

Félix d'Hérelle Reference Center for Bacterial Viruses

Last modifications: 2011-02-24