

Media recipe

CAY

Détails

For 1 liter of CAY:

1.5 g of casamino acids.

2.5 g of yeast extract.

20 ml of glycerol.

1 g of NaCl.

1 g of K₂HPO₄.

0.1 g of MgSO₄·7H₂O.

0.005 g of CaCO₃.

Complete to 1 liter with water.

For CAY-agar:

Add 15 g of agar per liter.

For CAY-top agar:

Add 7.5 g of agar per liter.

Autoclave.

Last modifications: 2011-01-24