

Media recipe

GAM

Détails

For 1 liter of GAM:

10 g of peptone.

3 g of soybean peptone.

10 g of proteose peptone.

13 g of digest serum powder.

5 g of yeast extract.

2,2 g of meat extract.

1,2 g of liver extract.

3 g of glucose.

2.5 g of dipotassium phosphate.

3 g of NaCl.

5 g of soluble starch.

0.3 g of L-cysteine-HCl.

0.3 g of sodium thioglycolate.

Complete to 1 liter with water.

For GAM-agar:

Add 12g of agar per liter.

For GAM-top agar:

Add 7.5 g of agar per liter.

Autoclave.

Last modifications: 2011-01-24